


I.O. DINE™ COMPLEX BOLUS

Water Soluble Bolus for Flushing and Topical Disinfection

When you need one product to flush, disinfect and cleanse, I.O. Dine Complex Bolus is your answer. This Bolus is a practical, flexible and cost-effective antiseptic that suits many of your veterinary practice needs.

- Polyvinylpyrrolidone – Iodine Bolus
- Buffered urea base
- When dissolved in water can be use as a:
 - Topical disinfectant
 - External wound flush
 - General flush
- Essentially non-irritating
- Each bolus contains 250 mg of Polyvinylpyrrolidone



Why I.O. Dine Bolus?

- Formulated for a variety of uses
- Quick dissolving for easy mixing
- Suitable for use with all animals



I.O. DINE™ COMPLEX BOLUS

GUARANTEED ANALYSIS:

250 mg. polyvinylpyrrolidone iodine complex buffered urea.

CONTAINS

250 mg. available (active) iodine as polyvinylpyrrolidone iodine complex in a buffered urea base.

MIXING INSTRUCTIONS

Dissolve one bolus in 50 to 100 mL of water. Apply solution as a soak, flush or saturate gauze with mixture and cleanse affected area.

STORAGE AND SAFETY

Store in a cool dry area

Keep out of reach of children

Use only as directed

Not for human use

For veterinary use only.

